### <u>Unity</u> <u>Seahaven</u>

| Script                 | Important Variables  |
|------------------------|--|
| VRTK_Touchpad walking  | <ul><li>Footstep sounds</li><li>Walking speed</li></ul>  |
| VRTK_headset Collision | • Ignore Target with Tag X   |
| Pupil Gaze Tracker     | Pupilrecording     EyesOpen     RayDistance     Server IP     Service Port     trainingStarted |
| PupilCalibMarker       | Points for validation  |
| EyeGazeRenderer        | •_image.enabled  |
| recorder               | IsRec     VPNum (counts automatically, only specify when needed)                               |
| Screenshot             | My Camera     ResWidthN     resHightN  |
| Auto Intensity         | Day rotate speed     All other parameters you want to play with                                |
| (ShowMap)              | Not in Seahven 2.0   |

### Commands: To be used during a session

| Effect  | Associated Script(s)  |
|---|---|
| Start/Stop pupil data recording   | PupilGazeTracker  |
| Start Calibration (17 points, whole field)  | PupilCalibMarker -> PupilGazeTracker  |
| Start 2D Validation (9 points, central)   | PupilCalibMarker -> PupilGazeTracker  |
| Start 3D Validation (9 points, central)   | PupilCalibMarker -> PupilGazeTracker  |
| Stop Callibration & Validation  | PupilCalibMarker -> PupilGazeTracker  |
| Stop all recordings, save them and quit game  | PupilGazeTracker, recorder  |
| Start VR Training -> stat pupil recording + Imoty recording, transform position to start position | PupilGazeTracker -> recorder  |
| Fast validation with one point  | PupilGazeTracker -> recorder  |
| Pause the session and all recordings  | PupilGazeTracker -> recorder  |
|   | Start/Stop pupil data recording Start Calibration (17 points, whole field) Start 2D Validation (9 points, central) Start 3D Validation (9 points, central) Stop Callibration & Validation Stop all recordings, save them and quit game Start VR Training -> stat pupil recording + Imoty recording, transform position to start position Fast validation with one point |

### DrawViewingPath:

| Script            | Important Variables  |
|-------------------|--|
| DrawViewingPath * | Randomize     VPNum     RayDistance  |
| Draw All Paths    | Generates all ViewedHouses Scripts of the *NumAnalyze* first recordings in your folder subjetcs in your recording folder (Don't do too many at once, it takes a while) |

# Recordings:

| File Name                                   | Data   | Script           |
|---|--|------------------|
| EyesOnScreen_VP#                            | Variable: Gazes = 2D coordinates of gaze (normalized) =(CenterX,CenterY) or (0.000000, 0.000000) | PupilGazeTracker |
| EyeBoxPos_VP#                               | Variable: BoxPos = 3D coordinates of box position  | PupilGazeTracker |
| Positions_VP#                               | (x,y,z,rx,ry,rz,timestamp (in sec),PupilTimeStamp)   | Recorder         |
| Validation2D+Num                            | Degree of error for each point + avg + time + last cal + error in x and y dir (+avg)             | PupilGazeTracker |
| Validation3D+Num                            | Degree of error for each point + avg + time + last cal   | PupilGazeTracker |
| (MapViews)                                  | Duration of each time the map was looked at  | ShowMap          |
| ViewedHouses_VP#                            | HouseViewed, distance, timestamp (sec. since start)  | DrawViewingPath  |
| Heatmap3D(Random<br>Gaze/RandomPos)_V<br>P# | List of (x,y,distance) of fixations during the session   | DrawViewingPath  |

The top 6 recordings are being automatically made during a session (MapViews are currently not used since there is no minimap in the experimental setup). To create Viewed houses and heatmap data open the Unity project 'DrawViewingPath', put in the VP number of the subject of interest into the respective field in the inspector and press play. After some calculation time you will be able to walk through the 3D city model and see all Gaze hit points. The two.txt files with the recordings have been created automatically. Additionally there are two checkboxes where you can pick a randomization (position or gaze randomization).

## • ViewedHouses:

House# -> House was looked at NH -> No house was looked at Distance =

0 -> Low confidence (<0.5) 200 -> No object hit (eg: sky) D -> House/obj. hit

## Matlab: (for new data format of Seahaven 2.0)

All important variables can be set on the top of each sript. Outputs are automatically saved. Make sure to have all the files named in the column 'Input' in the folder which you access with the respective MATLAB script. If you don't have a certain input file you might have to run DrawViewingPath (explained above) or a nother MATLAB script first. Check the 'Output' column to see which output is created by which script.

| Script             | Input   | Variables                           | Output   |
|--------------------|---|-------------------------------------|--|
| ValidationAnalysis | All Validation_VP#_Val#.txt files                                 | Path: Where should results be saved | validations: struct with all validation info for all subjects     Prints some overall analysis                         |
| PositionAnalysis   | Positions_VP#.txt   | • sourcepath                        | Map_VP_#.mat     North_VP_#.mat     Path_VP_#.mat  |
| Analysis_Map       | <ul><li>Path_VP_#</li><li>North_VP_#</li><li>(map_VP_#)</li></ul> | • savepath                          | OverlaidMap.jpeg IndividualNorth.jpeg SubjectHeatmap.jpeg (save images manually) Repeated measurements Path comparison |
| Heatmap3D          | 3DHeatmap(RandomX)_VP#.txt  | • VPNum                             | Heatmap (.jpeg)  |

|                       |  | • Condition<br>• Savepath                                       | .mat file of x,y,d,c (c=density at point)  |
|-----------------------|--|---|--|
| AnalyzeAllViews       | ViewedHouses_VP                                  | • savepath  | TimeLine (.jpeg) NumViewsD (.mat)  |
| Analysis_ViewedHouses | NumViewsD.mat                                    | • PartList<br>• savepath  | TotalNum{VPRange}.mat ViewingStats{VPRange}.mat  |
| PerformanceAnalysis   | ViewedHouses_VP#.txt<br>AlignmentVR_SubjNo_#.mat | Dname: directory of viewed houses                               | - Plot for Performance vs. viewing time/distance/d.variance                                      |
| GazeStandVSWalk       | • EyesOnScreen_VP#.txt • Positions_VP#.txt       | PartList     sourcepath   | GazeWalkStand{VPRange} (.jpeg)     Variances{VPRange} (.mat)                                     |
| LeftRightTurns        | EyesOnScreen_VP#.txt     Positions_VP#.txt       | PartList     Sourcepath     IntervalLength     TurnSignificance | GazeLeftRight{VPRange} (.jpeg)     TtestsLR{VPRange}.mat   |
| Entropy               | ViewedHouses_VP#     AlignmentVR_SubjNo_#.mat    | PartList     SourcPath     Intervallength                       | Entropy_IntervalLen_NumSJs(.mat) -> all entropies     Performance_EntropyXsecIntervalYSJs(.jpeg) |